




CPAP MASK LINERS

Company	ATHFAR 828 Inc	PAD A CHEEK LLC	WakeWell
Mask Liner	 Airo CPAP Mask Liner	 Mask Liner Full Face D	 Snugz Full Face Mask Liners
Website	Not provided	www.padacheek.com	www.sleepsnugz.com/buy-it-now/full-face
Description	Designed for hospitals and sleep labs to assist with many of the problems experienced by CPAP users. Fast for medical professionals to affix and remove; creates a better seal between mask and face; reduces noisy air leaks; minimizes skin irritations; absorbs moisture.	Each of our liners is designed for a specific mask. Soft, smooth, washable, durable, and reusable, our liners secure around a mask and stay with it when it is removed. Many have a Nose Pad built into the design to protect the bridge of the nose. If you don't see what you need, PAD A CHEEK can design it.	Snugz Full Face Mask Liners are machine washable, one-size-fits-most liners designed for most mouth-and-nose masks. They help reduce noisy leaks, prevent redness and irritation, and feel great against your skin.
Compatible CPAP Masks	full face masks	ResMed AirFit F20/AirTouch F20	All sizes of all nose-and-mouth masks except for the ResMed Mirage Quattro, the extra-large Fisher and Paykel Forma, and the large Resprionics ComfortGel. They are also not recommended for under-the-nose masks such as the Resprionics Amara View or DreamWear and the ResMed F30, or total face masks such as the Resprionics FitLife.
Materials	98% cotton, 2% Lycra	microfiber	nylon, polyester, Lycra
Fitting & Attachment	Available in Small, Medium, and Large. Three tabs affix directly to the mask.	Remove the headgear, line up the teardrop-shaped hole over the nose/mouth hole of the mask. Slide the arms of the mask into the slots of the liner. Secure with the Velcro.	Stretch the liner to fit over the mask.
Cost	Not provided	\$23	\$30 for 2
Cleaning	N/A (disposable)	Hand wash in hot soapy water; rinse well. Can also be sanitized in a UV or ozone cleaner.	machine washable
Estimated Useful Life	12 hours	6 months to a year or longer	45 days (including machine washes between uses)
Sample Request (for Sleep Professionals, not Patients)	Contact customer-care@athfar828.com	Contact karen@padacheek.com	Contact info@sleepsnugz.com
Become a Reseller (for Sleep Facilities)	N/A	Contact karen@padacheek.com for a reseller form	Contact info@sleepsnugz.com ●

Information based on data submitted by CPAP mask liner marketers. Sleep Review strives for accuracy in all data but cannot be held responsible for claims made by marketers. All CPAP mask liners may not be included. Email editor@sleepreviewmag.com to be considered for the next update.